



RPL3L Antibody

Product Code	CSB-PA040234
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92901
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein L3L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPL3L; 60S ribosomal protein L3-like
Immunogen Species	Homo sapiens (Human)
Target Names	RPL3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.