





## **RPL10** Antibody

| Product Code               | CSB-PA040238                                                                                                          |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Storage                    | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                         |
| Uniprot No.                | P27635                                                                                                                |
| Immunogen                  | Synthesized peptide derived from the Internal region of Human QM.                                                     |
| Raised In                  | Rabbit                                                                                                                |
| Species Reactivity         | Human, Mouse, Rat                                                                                                     |
| <b>Tested Applications</b> | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000                                                                             |
| Form                       | Liquid                                                                                                                |
| Conjugate                  | Non-conjugated                                                                                                        |
| Storage Buffer             | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Purification Method        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype                    | IgG                                                                                                                   |
| Alias                      | RPL10; DXS648E; QM; 60S ribosomal protein L10; Laminin receptor homolog; Protein QM; Tumor suppressor QM              |
| Immunogen Species          | Homo sapiens (Human)                                                                                                  |
| Target Names               | RPL10                                                                                                                 |