





RPL22 Antibody

Product Code	CSB-PA040241
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35268
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein L22.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
	The antibody was affinity-purified from rabbit antiserum by affinity-
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Purification Method Isotype	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. IgG RPL22; 60S ribosomal protein L22; EBER-associated protein; EAP; Epstein-