









## RPL37A Antibody

<b>Product Code</b>	CSB-PA040244
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61513
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein L37A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
<b>Tested Applications</b>	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPL37A; 60S ribosomal protein L37a
Immunogen Species	Homo sapiens (Human)
Target Names	RPL37A