





RPLP2 Antibody

Product Code	CSB-PA040245
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05387
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein LP2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
	ememategraph, demig spriepe speems immunegem
Isotype	IgG
Isotype Alias	
	IgG RPLP2; D11S2243E; RPP2; 60S acidic ribosomal protein P2; Renal carcinoma