



RPS5 Antibody

Product Code	CSB-PA040246
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46782
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein S5.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS5; 40S ribosomal protein S5
Immunogen Species	Homo sapiens (Human)
Target Names	RPS5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.