



Phospho-RPS6 (S235/S236) Antibody

Product Code	CSB-PA040248
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62753
Immunogen	Synthesized peptide derived from Human Ribosomal Protein S6 around the phosphorylation site of S235/S236.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS6; OK/SW-cl.2; 40S ribosomal protein S6; Phosphoprotein NP33
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.