




Phospho-RPS6KA1 (S380) Antibody

Product Code	CSB-PA040251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15418
Immunogen	Synthesized peptide derived from Human Rsk-1 around the phosphorylation site of S380.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6 kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K; MAP kinase-activated protein kinase 1a; MAPK-activated protein kinase 1a; MAPKAP kinase 1a; MAPKAP
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA1

Image



Western Blot analysis of 293 cells using Phospho-Rsk-1 (S380) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.