



Phospho- RPS6KA1/RPS6KA3/RPS6KA2/RPS6KA6 (S221/227/S218/232) Antibody

Product Code	CSB-PA040255
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15418/P51812/Q15349/Q9UK32
Immunogen	Synthesized peptide derived from Human Rsk-1/2/3/4 around the phosphorylation site of S221/227/S218/232.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6 kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K; MAP kinase-activated protein kinase 1a; MAPK-activated protein kinase 1a; MAPKAP kinase 1a; MAPKAP
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA1/RPS6KA3/RPS6KA2/RPS6KA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.