



RPS20 Antibody

Product Code	CSB-PA040274
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60866
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein S20.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS20; 40S ribosomal protein S20
Immunogen Species	Homo sapiens (Human)
Target Names	RPS20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.