



S100A11 Antibody

Product Code	CSB-PA040286
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31949
Immunogen	Synthesized peptide derived from the Internal region of Human Calcizzarin.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	S100A11; MLN70; S100C; Protein S100-A11; Calcizzarin; Metastatic lymph node gene 70 protein; MLN 70; Protein S100-C; S100 calcium-binding protein A11
Immunogen Species	Homo sapiens (Human)
Target Names	S100A11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.