




SDC4 Antibody

Product Code	CSB-PA040308
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31431
Immunogen	Synthesized peptide derived from Human Syndecan-4 around the non-phosphorylation site of S179.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SDC4; Syndecan-4; SYND4; Amphiglycan; Ryudocan core protein
Immunogen Species	Homo sapiens (Human)
Target Names	SDC4

Image



Western Blot analysis of HepG2 cells using Syndecan-4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.