



# MRPL17 Antibody

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA050001                                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                         |
| <b>Uniprot No.</b>         | Q9NRX2                                                                                                                |
| <b>Immunogen</b>           | Synthesized peptide derived from the C-terminal region of Human MRP-L17.                                              |
| <b>Raised In</b>           | Rabbit                                                                                                                |
| <b>Species Reactivity</b>  | Human                                                                                                                 |
| <b>Tested Applications</b> | IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000                                                                            |
| <b>Form</b>                | Liquid                                                                                                                |
| <b>Conjugate</b>           | Non-conjugated                                                                                                        |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Alias</b>               | MRPL17; LIP2; 39S ribosomal protein L17; mitochondrial; L17mt; MRP-L17; LYST-interacting protein 2                    |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                  |
| <b>Target Names</b>        | MRPL17                                                                                                                |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                           |