



BLM Antibody

Product Code	CSB-PA050013
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54132
Immunogen	Synthesized peptide derived from Human BLM around the non-phosphorylation site of T99.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BLM; RECQ2; RECQL3; Bloom syndrome protein; DNA helicase; RecQ-like type 2; RecQ2; RecQ protein-like 3
Immunogen Species	Homo sapiens (Human)
Target Names	BLM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.