



SRSF5 Antibody

Product Code	CSB-PA050022
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13243
Immunogen	Synthesized peptide derived from the Internal region of Human SRp40.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SRSF5; HRS; SFRS5; SRP40; Serine/arginine-rich splicing factor 5; Delayed-early protein HRS; Pre-mRNA-splicing factor SRP40; Splicing factor; arginine/serine-rich 5
Immunogen Species	Homo sapiens (Human)
Target Names	SRSF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.