




Phospho-SGK1 (S422) Antibody

Product Code	CSB-PA050044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00141
Immunogen	Synthesized peptide derived from Human SGK1 around the phosphorylation site of S422.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SGK1; SGK; Serine/threonine-protein kinase Sgk1; Serum/glucocorticoid-regulated kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	SGK1

Image



Western Blot analysis of VEC HY926 cells using Phospho-SGK1 (S422) Polyclonal Antibody