



SGK1 Antibody

Product Code	CSB-PA050046
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00141
Immunogen	Synthesized peptide derived from Human SGK1 around the non-phosphorylation site of T256.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SGK1; SGK; Serine/threonine-protein kinase Sgk1; Serum/glucocorticoid-regulated kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	SGK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.