



ANAPC1 Antibody

Product Code	CSB-PA050059
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1A4
Immunogen	Synthesized peptide derived from Human APC1 around the non-phosphorylation site of S355.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANAPC1; TSG24; Anaphase-promoting complex subunit 1; APC1; Cyclosome subunit 1; Mitotic checkpoint regulator; Testis-specific gene 24 protein
Immunogen Species	Homo sapiens (Human)
Target Names	ANAPC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.