



# ANAPC1 Antibody

|                            |                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA050059                                                                                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                 |
| <b>Uniprot No.</b>         | Q9H1A4                                                                                                                                        |
| <b>Immunogen</b>           | Synthesized peptide derived from Human APC1 around the non-phosphorylation site of S355.                                                      |
| <b>Raised In</b>           | Rabbit                                                                                                                                        |
| <b>Species Reactivity</b>  | Human, Mouse                                                                                                                                  |
| <b>Tested Applications</b> | IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000                                                                                                    |
| <b>Form</b>                | Liquid                                                                                                                                        |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                       |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                         |
| <b>Isotype</b>             | IgG                                                                                                                                           |
| <b>Alias</b>               | ANAPC1; TSG24; Anaphase-promoting complex subunit 1; APC1; Cyclosome subunit 1; Mitotic checkpoint regulator; Testis-specific gene 24 protein |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                          |
| <b>Target Names</b>        | ANAPC1                                                                                                                                        |