





Phospho-NCF1 (S345) Antibody

Product Code	CSB-PA050086
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14598
Immunogen	Synthesized peptide derived from Human p47-phox around the phosphorylation site of S345.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NCF1; NOXO2; SH3PXD1A; Neutrophil cytosol factor 1; NCF-1; 47 kDa autosomal chronic granulomatous disease protein; 47 kDa neutrophil oxidase factor; NCF-47K; Neutrophil NADPH oxidase factor 1; Nox organizer 2; Nox-organizing protein 2; SH3
Immunogen Species	Homo sapiens (Human)
Target Names	NCF1

Image



Western Blot analysis of Jurkat cells using Phospho-p47-phox (S345) Polyclonal Antibody



Western blot analysis of 3T3 lysis using Phospho-p47-phox (S345) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.