



FOXD4/FOXD4L4/FOXD4L5/FOXD4L2/FOXD4L3/ FOXD4L6/ Antibody

Product Code	CSB-PA050093
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12950/Q8WXT5/Q5VV16/Q6VB85/Q6VB84/Q3SYB3
Immunogen	Synthesized peptide derived from the Internal region of Human FoxD4/D4L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXD4; FKHL9; FOXD4A; FREAC5; Forkhead box protein D4; Forkhead-related protein FKHL9; Forkhead-related transcription factor 5; FREAC-5; Myeloid factor-alpha; FOXD4L4; FOXD4B; Forkhead box protein D4-like 4; FOXD4-like 4; Forkhead box prote
Immunogen Species	Homo sapiens (Human)
Target Names	FOXD4/FOXD4L4/FOXD4L5/FOXD4L3/FOXD4L6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.