



BMP8B/BMP8A Antibody

Product Code	CSB-PA050097
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34820/Q7Z5Y6
Immunogen	Synthesized peptide derived from the Internal region of Human BMP-8.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BMP8B; BMP8; Bone morphogenetic protein 8B; BMP-8; BMP-8B; Osteogenic protein 2; OP-2; BMP8A; Bone morphogenetic protein 8A; BMP-8A
Immunogen Species	Homo sapiens (Human)
Target Names	BMP8B/BMP8A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.