




SLC15A1 Antibody

Product Code	CSB-PA050098
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46059
Immunogen	Synthesized peptide derived from the Internal region of Human PEPT1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SLC15A1; PEPT1; Solute carrier family 15 member 1; Intestinal H(+)/peptide cotransporter; Oligopeptide transporter; small intestine isoform; Peptide transporter 1
Immunogen Species	Homo sapiens (Human)
Target Names	SLC15A1

Image



Western Blot analysis of MOUSE-BRAIN
MOUSE-HEART MOUSE-LUNG MOUSE-
KIDNEY cells using PEPT1 Polyclonal Antibody