




Phospho-BMX (Y566) Antibody

Product Code	CSB-PA050105
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51813
Immunogen	Synthesized peptide derived from Human Bmx around the phosphorylation site of Y566.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BMX; Cytoplasmic tyrosine-protein kinase BMX; Bone marrow tyrosine kinase gene in chromosome X protein; Epithelial and endothelial tyrosine kinase; ETK; NTK38
Immunogen Species	Homo sapiens (Human)
Target Names	BMX

Image



Western blot analysis of 4T1 using p-Bmx (Y566) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.