




SP1 Antibody

Product Code	CSB-PA050124
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08047
Immunogen	Synthesized peptide derived from Human Sp1 around the non-phosphorylation site of T453.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SP1; TSFP1; Transcription factor Sp1
Immunogen Species	Homo sapiens (Human)
Target Names	SP1

Image



Western Blot analysis of A549 cells using Sp1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.