



Phospho-BRCA1 (S1524) Antibody

Product Code	CSB-PA050136
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38398
Immunogen	Synthesized peptide derived from Human BRCA1 around the phosphorylation site of S1524.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53
Immunogen Species	Homo sapiens (Human)
Target Names	BRCA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.