




BRCA1 Antibody

Product Code	CSB-PA050138
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38398
Immunogen	Synthesized peptide derived from Human BRCA1 around the non-phosphorylation site of S1457.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53
Immunogen Species	Homo sapiens (Human)
Target Names	BRCA1

Image



Western Blot analysis of 293 cells using BRCA1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.