



Phospho-SRF (S103) Antibody

Product Code	CSB-PA050141
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11831
Immunogen	Synthesized peptide derived from Human SRF around the phosphorylation site of S103.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SRF; Serum response factor; SRF
Immunogen Species	Homo sapiens (Human)
Target Names	SRF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.