





Phospho-BRAF (S446) Antibody

Product Code	CSB-PA050148
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15056
Immunogen	Synthesized peptide derived from Human Raf-B around the phosphorylation site of S446.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BRAF; BRAF1; RAFB1; Serine/threonine-protein kinase B-raf; Proto-oncogene B-Raf; p94; v-Raf murine sarcoma viral oncogene homolog B1
Immunogen Species	Homo sapiens (Human)
Target Names	BRAF

Image



Western Blot analysis of 293 cells using Phospho-Raf-B (S446) Polyclonal Antibody



Western blot analysis of K562 293T lysis using Phospho-Raf-B (S446) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.