



Phospho-BRCA2 (S3291) Antibody

Product Code	CSB-PA050157
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51587
Immunogen	Synthesized peptide derived from Human BRCA2 around the phosphorylation site of S3291.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BRCA2; FACD; FANCD1; Breast cancer type 2 susceptibility protein; Fanconi anemia group D1 protein
Immunogen Species	Homo sapiens (Human)
Target Names	BRCA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.