



Phospho-STAT2 (Y631) Antibody

Product Code	CSB-PA050163
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52630
Immunogen	Synthesized peptide derived from Human Stat2 around the phosphorylation site of Y631.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAT2; Signal transducer and activator of transcription 2; p113
Immunogen Species	Homo sapiens (Human)
Target Names	STAT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.