




STAT5A/STAT5B Antibody

Product Code	CSB-PA050170
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42229/P51692
Immunogen	Synthesized peptide derived from Human Stat5a/b around the non-phosphorylation site of S726/731.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)

Image



Western blot analysis of 293T lysis using Stat5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.