




Phospho-STK11 (T189) Antibody

Product Code	CSB-PA050177
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15831
Immunogen	Synthesized peptide derived from Human LKB1 around the phosphorylation site of T189.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STK11; LKB1; PJS; Serine/threonine-protein kinase STK11; Liver kinase B1; LKB1; hLKB1; Renal carcinoma antigen NY-REN-19
Immunogen Species	Homo sapiens (Human)
Target Names	STK11

Image



Western Blot analysis of 3T3 cells using Phospho-LKB1 (T189) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.