



Phospho-SYK (Y352) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA050187 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P43405 |
| Immunogen | Synthesized peptide derived from Human Syk around the phosphorylation site of Y352. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SYK |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |