



Phospho-SYK (Y323) Antibody

Product Code	CSB-PA050188
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43405
Immunogen	Synthesized peptide derived from Human Syk around the phosphorylation site of Y323.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk
Immunogen Species	Homo sapiens (Human)
Target Names	SYK