




TACR1 Antibody

Product Code	CSB-PA050201
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25103
Immunogen	Synthesized peptide derived from the Internal region of Human NK-1R.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TACR1; NK1R; TAC1R; Substance-P receptor; SPR; NK-1 receptor; NK-1R; Tachykinin receptor 1
Immunogen Species	Homo sapiens (Human)
Target Names	TACR1

Image



Western Blot analysis of MOUSE-BRAIN COLO cells using NK-1R Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.