




Phospho-MAP3K7 (T187) Antibody

Product Code	CSB-PA050206
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43318
Immunogen	Synthesized peptide derived from Human Tak1 around the phosphorylation site of T187.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP3K7; TAK1; Mitogen-activated protein kinase kinase kinase 7; Transforming growth factor-beta-activated kinase 1; TGF-beta-activated kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K7

Image



Western Blot analysis of 293T cells using Phospho-Tak1 (T187) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.