





TBX3 Antibody

Product Code	CSB-PA050217
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15119
Immunogen	Synthesized peptide derived from the Internal region of Human TBX3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TBX3; T-box transcription factor TBX3; T-box protein 3
Immunogen Species	Homo sapiens (Human)
Target Names	TBX3

Image



Western Blot analysis of KB cells using TBX3 Polyclonal Antibody



Western blot analysis of mouse-kidney mouse-brain KB 293T lysis using TBX3 antibody.

