




VPS72 Antibody

Product Code	CSB-PA050222
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15906
Immunogen	Synthesized peptide derived from the Internal region of Human VPS72.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VPS72; TCFL1; YL1; Vacuolar protein sorting-associated protein 72 homolog; Protein YL-1; Transcription factor-like 1
Immunogen Species	Homo sapiens (Human)
Target Names	VPS72

Image



Western blot analysis of 293T lysis using VPS72 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.