




TERF1 Antibody

Product Code	CSB-PA050240
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54274
Immunogen	Synthesized peptide derived from Human TRF1 around the non-phosphorylation site of S219.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TERF1; PIN2; TRBF1; TRF; TRF1; Telomeric repeat-binding factor 1; NIMA-interacting protein 2; TTAGGG repeat-binding factor 1; Telomeric protein Pin2/TRF1
Immunogen Species	Homo sapiens (Human)
Target Names	TERF1

Image



Western Blot analysis of 293 cells using TRF1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.