



Phospho-TGFBR2 (S225) Antibody

Product Code	CSB-PA050249
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37173
Immunogen	Synthesized peptide derived from Human TGFβ RII around the phosphorylation site of S225.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TGFBR2; TGF-beta receptor type-2; TGFR-2; TGF-beta type II receptor; Transforming growth factor-beta receptor type II; TGF-beta receptor type II; TbetaR-II
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.