



TGFBR3 Antibody

Product Code	CSB-PA050252
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03167
Immunogen	Synthesized peptide derived from the Internal region of Human TGFβ RIII.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TGFBR3; Transforming growth factor beta receptor type 3; TGF-beta receptor type 3; TGFR-3; Betaglycan; Transforming growth factor beta receptor III; TGF-beta receptor type III
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.