




THBS4 Antibody

Product Code	CSB-PA050258
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35443
Immunogen	Synthesized peptide derived from the Internal region of Human Thrombospondin 4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	THBS4; TSP4; Thrombospondin-4
Immunogen Species	Homo sapiens (Human)
Target Names	THBS4

Image



Western blot analysis of 293T lysis using Thrombospondin 4 antibody.