



Phospho-TK1 (S13) Antibody

Product Code	CSB-PA050266
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04183
Immunogen	Synthesized peptide derived from Human Thymidine Kinase around the phosphorylation site of S13.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TK1; Thymidine kinase; cytosolic
Immunogen Species	Homo sapiens (Human)
Target Names	TK1