




Phospho-TOP2A (S1525) Antibody

Product Code	CSB-PA050286
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11388
Immunogen	Synthesized peptide derived from Human Topo II α around the phosphorylation site of S1525.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TOP2A; TOP2; DNA topoisomerase 2-alpha; DNA topoisomerase II; alpha isozyme
Immunogen Species	Homo sapiens (Human)
Target Names	TOP2A

Image



Western blot analysis of 293T 3T3 COLO using p-Topo II α (S1525) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.