



Phospho-TP53BP1 (S25) Antibody

Product Code	CSB-PA050308
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12888
Immunogen	Synthesized peptide derived from Human 53BP1 around the phosphorylation site of S25.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1
Immunogen Species	Homo sapiens (Human)
Target Names	TP53BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.