



TP53BP1 Antibody

Product Code	CSB-PA060001
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12888
Immunogen	Synthesized peptide derived from Human 53BP1 around the non-phosphorylation site of S1778.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1
Immunogen Species	Homo sapiens (Human)
Target Names	TP53BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.