



UBE2G1 Antibody

Product Code	CSB-PA060036
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62253
Immunogen	Synthesized peptide derived from the Internal region of Human UBE2G1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	UBE2G1; UBE2G; Ubiquitin-conjugating enzyme E2 G1; E217K; UBC7; Ubiquitin carrier protein G1; Ubiquitin-protein ligase G1
Immunogen Species	Homo sapiens (Human)
Target Names	UBE2G1