



Phospho-VASP (S238) Antibody

Product Code	CSB-PA060045
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50552
Immunogen	Synthesized peptide derived from Human VASP around the phosphorylation site of S238.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VASP; Vasodilator-stimulated phosphoprotein; VASP
Immunogen Species	Homo sapiens (Human)
Target Names	VASP