



Phospho-VAV1 (Y174) Antibody

Product Code	CSB-PA060053
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15498
Immunogen	Synthesized peptide derived from Human Vav1 around the phosphorylation site of Y174.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VAV1; VAV; Proto-oncogene vav
Immunogen Species	Homo sapiens (Human)
Target Names	VAV1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.